Therabody

# RecoveryTherm

Hot & Cold Vibration Knee

Powered by *Cryothermal* Technology

User Manual

## Table of contents

EN	3-10
ES	11-18
FR	19-26
IT	27-34
DE	35-42
FRCA	43-50
SC	51-58

## Product Overview

# RecoveryTherm

EN

## Hot & Cold Vibration Knee

Precise cold, heat, and vibration for advanced knee therapy.

From the creators of Theragun comes RecoveryTherm Hot & Cold Vibration Knee, precisely delivering heating, cooling, and vibration therapy in one wearable garment, the first of its kind. From improving athletic performance to enhancing muscle recovery and general wellness, the RecoveryTherm Hot & Cold Vibration Knee offers a solution for everyone.

### Intended Use

The RecoveryTherm Hot & Cold Vibration Knee is intended to help alleviate pain, improve athletic performance, and enhance muscle recovery post-workout. While cooling, heating, and vibration offer unique therapeutic benefits, rapidly alternating between the three, known as contrast therapy, optimizes and increases blood flow to reduce pain and soreness.

# Getting started

## EN Step 1

### Put on your garment

1. Find a comfortable position sitting down, with your legs either at a 90-degree angle or fully extended.
2. Unbuckle the straps and extend the garment.
3. To put the RecoveryTherm Hot & Cold Vibration Knee on, place the garment's center open circle on your knee cap. The blue fans should be facing outwards, and the control panel should be resting on your thigh. Ensure that nothing is blocking the airflow to the blue fans.
4. Secure the RecoveryTherm Hot & Cold Vibration Knee by buckling each strap. The straps should provide compression without being uncomfortable or painful. If standing up, readjust the top strap for increased compression. Add the extenders to the end of the straps when using the device for larger muscles.
5. For additional support, attach the third strap (found inside the pouch in the carrying case). To attach, Velcro one of the sides using the rubber connectors that are to the right and left of the knee cap. Use the third strap in instances where further support is needed, for example, when standing or moving during treatment. The use of the third strap is optional and fully removable.
6. For safe and optimal use, store your garment in a cool, dry location (below 78 degrees Fahrenheit/26 degrees Celsius) to protect the temperature of the metal plates.

Scan the QR code to get started.

## Step 2

### Use your Garment

1. Insert the battery into the control panel. The Therabody logo should be facing outwards, with the charging port on the bottom. Battery insertion should be seamless and not forced.
2. Press and hold the power button for two seconds.
3. Select one of three treatment options (cold, hot, or contrast) using the respective mode buttons.
4. Short press the snowflake button to use the cold treatment only. The blue lights around the Therabody logo will illuminate when the cold setting is on. The vibration will automatically start on low; short press the vibration button to toggle between patterns (low, high, wave, off). The treatment will run for approximately 20 minutes. The LEDs will turn off sequentially as each treatment mode is complete; e.g., the blue lights will turn off after a cooling mode. Once the treatment is complete, the device will turn off.
5. Short press the fire button to use the hot treatment only. The orange lights around the Therabody logo will illuminate when the hot setting is on. The vibration will automatically start on high; short press the vibration button to toggle between patterns (low, high, wave, off). The treatment will run for approximately 15 minutes. Once the treatment is complete, the device will turn off.
6. Short press the contrast button to begin contrast therapy treatment. Both the orange and blue lights around the Therabody logo will illuminate when the contrast setting is on; blue indicates cold therapy, and orange indicates hot therapy. The vibration will automatically start on the wave setting; short press the vibration button to toggle between patterns (low, high, wave, off). The device will move from cold to hot for approximately 18 minutes. Once the treatment is complete, the device will turn off.

7. Vibration therapy can be used without hot, cold, or contrast therapy. Short press the vibration button once the device is turned on. Vibration can be added or customized at any point during your treatment; short press the vibration button to toggle between patterns (low intensity, high intensity, and wave).
8. Depending on the selected mode, treatment times can vary between 12-18 minutes.
9. To access the Pro Mode, short press the snowflake button three times. The Pro Mode is a cold treatment (with lower temperatures, reaching 8 degrees C) and a low vibration therapy.

### Step 3

#### Feel the difference. Every day.

Whether used for targeted pain relief, post-workout muscle recovery, or general soreness, RecoveryTherm Hot & Cold Vibration Knee will help keep you feeling your best.

#### See Warnings for important safety instructions

Inner Fabric: 89% germanium-infused nylon, 11% spandex  
Outer Fabric: 84% nylon, 16% spandex

SPOT CLEAN ONLY  
SEE USER MANUAL FOR CLEANING INSTRUCTIONS  
DO NOT MACHINE WASH

MADE IN CHINA

# Getting to know your wearable device

EN


Figure 1

- A** *Wrap garment*
- B** *Removable battery*
- C** *USB-C Port*
- D** *Three adjustable straps (the third strap is optional and can be found inside the pouch)*
- E** *Strap extenders*

# Protocols Overview

EN

Table 1. Preset protocols (see Figure 2)

	Step 1	Step 2	Step 3	Step 4	Step 5	Step 6	Total
<b>Heating</b> Temperature gradually increases (Orange light around Therabody logo)	2 Minutes	2 Minutes	2 Minutes	2 Minutes	2 Minutes	2 Minutes	12 Minutes
<b>Cooling</b> Maintains cold temperature (Blue light around Therabody logo)	3 Minutes	3 Minutes	3 Minutes	3 Minutes	3 Minutes	3 Minutes	18 Minutes
<b>Contrast</b> Temperature alternates between heating and cooling (Orange and blue lights around Therabody logo)	3 minutes Cooling	3 minutes Cooling	2 minutes Heating	3 minutes Cooling	3 minutes Cooling	2 minutes Heating	16 Minutes


Figure 2

## Table 2. Treatment type selection

When to use	What to Choose	Intended benefit	Pro Tip
<b>Pre-activity</b>	Hot	Warm-up joints Improve range of motion	Add vibration therapy to enhance blood flow.
<b>Post-activity</b>	Cold	Reduce joint swelling Reduce muscle soreness	Add vibration therapy to improve comfort during the cooling treatment.
<b>Post-activity General soreness</b>	Contrast	Reduce joint swelling Reduce muscle soreness	Add vibration therapy to enhance blood flow and improve comfort during the contrast treatment.

## Charging

EN

- To charge the RecoveryTherm Hot & Cold Vibration Knee, remove the battery from the control panel and connect it to the included USB-C adapter.
- The RecoveryTherm Hot & Cold Vibration Knee is USB-C enabled.
- The battery level is displayed in the LED indicator on the battery pack, next to the USB-C connector:
  - **Orange: low charge, blue: mid charge, green: fully charged.**
  - **When charging, the LED will flash to indicate it is properly connected and will turn solid green when complete.**

Note: If using an alternative adapter, ensure that it is from a certified manufacturer and has not suffered any structural damage. Do not try to use the device while charging.

## Cleaning Instructions:

1. The RecoveryTherm Hot & Cold Vibration Knee is NOT waterproof or machine washable. Do not allow excess liquids or water to enter the product.
2. Remove the battery before cleaning.
3. To clean, wipe the inner surface of the strap with a dampened microfiber cloth. Wipe down the metallic parts and the flexible fingers.

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