

Clinical Study Details

Clinical Study:

- This was a 12-week monadic clinical study on the efficacy of the Theraface PRO device. The study was completed by 35 study subjects who met the inclusion and non-inclusion criteria required. Study subjects were instructed to use the product starter guide and follow a weekly protocol for 12 weeks.

Participants:

- Female and male subjects, aged between 25 and 61 years old (mean age: 43 years old)
- Presenting all skin types, Fitzpatrick Skin Types I-IV
- Hispanic, African American, Caucasian, and Asian representation

Clinical Study Weekly Protocol:

| Monday | Tuesday | Wednesday | Thursday | Friday | Saturday | Sunday |
|--------------------------------|--------------------------------------|--------------------------------|--------------------------------------|--------------------------------|-------------------------------|--------------------------------|
| 6 min Percussive and LED | 6 min Cleanse and Microcurrent | 6 min Percussive and LED | 6 min Cleanse and Microcurrent | 6 min Percussive and LED | 6 min Microcurrent Only | 45 Min The PRO Treatment |

- **MWF** - Blue LED, Red LED, and Percussive
- **Tues & Thurs** - Cleanse and Microcurrent
- **Sat** - Microcurrent
- **Sun** - Cleanse, Microcurrent, Red LED, Red+IR LED, Percussive, Hot, Cold

TheraFace PRO is Clinically Shown to*:

- Improve elasticity, firmness, and facial contour*
- Improve the appearance of fine lines and wrinkles
- Improve skin tone evenness
- Improve overall healthy skin appearance
- Lift the neck and jowl lines.
- Improve facial contour
- Smooth skin texture
- Improve skin radiance/luminosity/glow
- Improve the appearance of sagging skin
- Improve elasticity, firmness, and skin extensibility

Clinical Analysis*:

- 97% of participants showed a reduction in wrinkles around the eyes
- 97% of participants had an increase in firmness
- 97% of participants showed an improvement in lifting
- 94% of participants experienced an improvement in skin texture and smoothness
- 94% of participants showed an improvement in skin tone evenness
- 91% of participants had improvements in fine lines around the eyes
- 87% of participants showed an improvement in elasticity
- 80% of participants had a lifting effect of the jowel*

* Data from 12-week clinical study on efficacy of TheraFace PRO device. Individual results may vary.

Consumer Perception Results*:

- **94%** reported their skin looked overall healthier.
- **91%** noticed an improvement in the overall appearance of their skin.
- **89%** observed a decrease in wrinkles and/or unhealthy-looking skin.
- **89%** felt their skin was more rejuvenated, radiant and had a glow.
- **89%** noticed an improvement in facial contour including improved muscle tone and skin tightness.
- **83%** noticed that their skin was less puffy.
- **83%** observed their skin looked and/or felt smoother and tighter.
- **80%** said their skin looks and/or feels firmer.

*Based on a 35-subject consumer perception study after 12-weeks of use.

